



KALEIDON probiotic is a food supplement based on the probiotic *Lactobacillus rhamnosus* GG (ATCC 53103) that benefits the balance of the intestinal flora.

KALEIDON probiotic comes in packets containing 10 X 1g sachets for oral use with a neutral flavour and wild berry flavour.

DOSE AND ADMINISTRATION METHOD:

1 Sachet a day.

Open the sachet and swallow directly in your mouth. Best taken with meals.

WARNINGS:

Keep out of the reach of children under three years of age.

Do not exceed the recommended daily dose.

Food supplements must not be considered as substitutes for a balanced and varied diet and a healthy lifestyle.

Excessive consumption may have laxative effects.

STORAGE METHOD:

The product contains probiotics which are sensitive to heat. It is recommended to avoid exposure to heat sources or temperature fluctuations and to store the product away from the light in a dry place preferably at a temperature below 25°C. The minimum shelf life refers to the product correctly stored in its unopened package.

INGREDIENTS:

Wild berry flavour.

Bulking agents: erythritol, xylitol; natural flavouring, *Lactobacillus rhamnosus* GG probiotic (ATCC 53103), coating agent: mono- and diglyceride fatty acids, anti-caking agent: magnesium oxide; coating agent: silicon dioxide.

Neutral flavour. Bulking agents: erythritol, xylitol; *Lactobacillus rhamnosus* GG (ATCC 53103), coating agent: mono- and diglyceride fatty acids, anti-caking agent: magnesium oxide, coating agent: silicon dioxide.

Gluten-free and naturally lactose free.

Ingredient

Average content in live cells per dose (1 sachet)

<i>Lactobacillus rhamnosus</i> GG	5x10 ⁹ CFU
-----------------------------------	-----------------------

One sachet contains not less than 5 billion live cells of *Lactobacillus rhamnosus* GG (ATCC 53103). This value is guaranteed until the end of the shelf life of the correctly stored product.

Privacy info: privacy section of the website www.firma-fi.it

F.I.R.M.A. S.p.A. Via di Scandicci 35-50143 FLORENCE, company in the Menarini Group.

Manufactured and packaged in the plant of Bogbindervej 6 – DK 3390 Hundested - Denmark